

FASTQ format - sequence ID line

```
@D7MHBFN1:202:D1BUDACXX:4:1101:1340:1967 1:N:0:CATGCA  
NATCTTCGGATCACTTTGGTCAAATTGAAACGATACAGAGAAGATTGTAAGTAACAATATTTACCAAGGTTCGAGTCATACTAACTCGTTGTCCTATAGT  
+  
#1=DDFFFHHHHHJJJJJJJHIJJJJJJIIJJGIIIJJJJJJIIJJJJJIIFGIIIIJJJJJJIEHJIIHHGFFF@?ADFEDDEDCDDDBDDBCDDDDEC
```

- D7MHBFI - unique instrument name
- 202 - run ID
- D1BUDACXX - flowcell ID
- 4 - flowcell lane
- 1101 - tile number within lane
- 1340 - x-coordinate of cluster within tile
- 1967 - y-coordinate of cluster within tile
- 1 - member of pair (1 or 2). Older versions: /1 and /2
- Y/N - whether the read failed quality control (Y = bad)
- 0 - none of the control bits are on
- CATGCA - index sequence (barcode)